

# OXLIFE INDEPENDENCE<sup>®</sup>

24/7 Solution. Reliable. Durable.



# OXLIFE INDEPENDENCE<sup>®</sup>

24/7 Solution. Reliable. Durable.



## COMPLETE SOLUTION

- 24/7 Solution - No Need for Stationary
- Two Billing Set Up Codes (E1390 & E1392) on One Device
- Five-year Warranty
- Quick Set Up in Hospital, Clinic or Home



## SUPERIOR PRODUCT FEATURES

- Continuous Flow to 3LPM and Twelve Pulse Settings
- Energy Smart Technology<sup>®</sup> = Longest Battery Life in Class
- Polycarbonate Structure with Best-In-Class Components
- Easy Battery Swap While Unit is in Operation

*“Now I am out there living my life. The Oxlife Independence<sup>®</sup> has given me my freedom back. Now I can go anywhere.”*

- John P., Patient

*“Thank you to O2 Concepts<sup>®</sup> for caring and hiring employees who reflect your service orientation. You are filling a need to provide people with compromised breathing the assistance to lead independent lives.”*

- Susan S., Care Provider



**BETTER**

## **PATIENT OUTCOMES**

- 24/7 single solution - continuous flow and large pulse dosages to ensure proper saturation - no need for stationary concentrator
- Easy set up in hospital or clinic provides convenient and efficient discharge for patients and providers - eliminates home delivery
- Proactive solution to reduce instances of COPD readmittance (one day in the hospital is more expensive than a year of oxygen therapy)



## PRODUCT SPECIFICATIONS<sup>1</sup>

Dimensions	Height: 20.2" • Width: 10.9" • Depth: 9.45"
Weight <sup>2</sup>	Device: 16.7 lbs • Battery: 1.35 lbs/each
Oxygen Flow Rate	Pulse: 0.5 - 6.0 in 0.5 increments (8 - 96ML) Continuous: 0.5 - 3 LPM in 0.5 increments
Battery Duration <sup>3</sup>	5.75 hours @ 2 pulse at 20 bpm 2.5 hours @ 2LPM continuous
Battery Recharge Time <sup>4</sup>	2.5 hours
Oxygen Concentration	87% - 95% at sea level at all settings
Oxygen Output Pressure	5 +/- 1 psig (34 +/- 7 kPa)
AC Power	Universal Input: 100 - 240 VAC 50 - 60 Hz • 2.0A maximum
DC Power	12 - 15 VDC • 15 amps
Operating Temperature	50°F (10°C) to 104°F (40°C)
Altitude Operating Range	0 - 13,123 feet (0 - 4,000 meters)
Transportation/Storage Humidity	0 - 95% non-condensing
Warranty	5 years on unit • 1 year on accessories

(1) O2 Concepts® reserves the right to revise specifications at anytime without notice

(2) Production weight excludes cart kit

(3) Battery duration presented is for extended battery system

(4) Battery recharge time increases if unit is operating while charging

## PULSE FLOW SETTINGS

SETTING	DOSAGE (ML)
0.5	8
1.0	16
1.5	24
2.0	32
2.5	40
3.0	48
3.5	56
4.0	64
4.5	72
5.0	80
5.5	88
6.0	96



**O<sub>2</sub> Concepts®**

6303 Waterford Blvd., Suite 150  
Oklahoma City, OK 73118  
877.867.4008 / O2-concepts.com



© 2017 O2 Concepts, LLC; all rights reserved; O2 Concepts®, Oxlife Independence® and Energy Smart Technology®, including certain related logos, are registered trademarks of O2 Concepts, LLC and may not be reproduced without express written consent.

900-1006 Rev: CE-C