


KSOC-10L & KSOC-5L OXYGEN CONCENTRATORS





- AMERICAN PSA TECHNOLOGY
- FRENCH IMPORTED MOLUCRE SIEVE BEDS
- DURABLE OIL FREE COMPRESSOR
- 24 HOUR CONTINUOUS WORKFLOW
- SELF-DIAGNOSTIC SYSTEM WITH MULTIPLE ALARM FUNCTIONS
 - Power off alarm
 - Overload protection
 - High and Low flow pressure alarm
 - Temperature alarm
 - Error Code indication
 - Oxygen purity alarm
- NEBULIZER PORT
- PURITY: >93% ±3 (AT ALL FLOW SETTINGS)
- COMPLY WITH THE WORLD HEALTH ORGANISATION REQUIREMENTS
- 2 YEAR WARRANTY
- INCLUDES A STARTER PACK



Large LED Screen including hours used and hours in use.

USA Technology, Oil free compressor



Low noise levels, light weight and high temperature protection

ACCESSORIES INCLUDED	
Description	Qty (pcs)
Oxygen concentrator	1
Warning paper	1
Foam filter	1
Oxygen cannula	2
Spanner	1
Humidifier bottle	1

Model	Flow Range	Oxygen Purity	Input Voltage	Output Pressure	Noise	Packing Size	G.W./N.W.	Power
KSOC-5	1-5L	93%±3%	220V/110V	30-70kPa	≤45dB	39×31×59cm	22kg/20kg	350W
KSOC-10	1-10L	93%±3%	220V/110V	30-70kPa	≤52dB	50×42×73cm	37kg/32kg	450W