


ALICE PDx



Alice® PDx™ is a diagnostic device for the portable diagnosis of cardio-respiratory sleep disorders to be used by labs, office based doctors and HST.

HIGHLIGHTS:

- 1) For Ambulatory Sleep Study (Home sleep study)
- 2) For standard patients with no co-morbidities
- 3) Limited channels: **EEG (1-3), EOG (1-2), EMG (3), ECG (3), Respiratory parameters (Thermistor, nasal cannula, thoracic & abdominal), SPO2, Body position, Snore microphone**
- 4) Therapy device connectivity
- 5) Versatile, Superior performance, Easy!
- 6) 20-24 hr Battery
- 7) 1 GB SD Card

ALICE PDx

FEATURES & BENEFITS:

FEATURES	BENEFITS
Good study indicator	<ul style="list-style-type: none">- Reduce technologist time- Limit unnecessary transportation to and from the patient's home
Auto Start / Stop	<ul style="list-style-type: none">- Start and Stop the acquisition whatever the patient is doing (forget to switch on/Off , etc..)- Can program several following acquisitions (up to SD card capacity)
From simple polygraphy to Portable sleep scoring (with options)	<ul style="list-style-type: none">-Simpler to handle: same device, same software.-Cost effective
Versatile	Can choose between basic PG, ECG Option , EXG option
Easy	<ul style="list-style-type: none">-LCD interface for Patient-Confirm Electrode Connection-Robust Carrying Case-Patients' Guide and Video-Battery and SD card check