



LOTUS MEDIC AIR

SPECIALIZING IN SLEEP APNEA THERAPY & HOME OXYGEN

Micomme i series

Sleep Apnea solution C2 C5 B5



Preheat



Separable
design



HIGH QUALITY
AFFORDABILITY


i Series Specifications

MODULE	C2	C5	B5
MODE	CPAP	CPAP & AUTO CPAP	CPAP & AUTO Bi-LEVEL
PRESSURE ACCURACY	± 0.2 cm H ₂ O	± 0.2 cm H ₂ O	± 0.2 cm H ₂ O
PRESSURE RANGE	4-20cm H ₂ O	4-20cm H ₂ O	4-20cm H ₂ O
RAMP TIME	0 to 45 min (5-minute increments)		
COMF PRESSURE RELIEF	1-3 levels		
HUMIDIFICATION LEVEL	1-5 levels (113 to 185°F/23 to 85°C		
RISE TIME	---	---	1-4 levels
SPLIT NIGHT	---	---	YES
DATA STORAGE CAPACITY	8G USB disk		
WEIGHT	1.72 kg		
MEAN SOUND LEVEL	≤ 30 dB		


Safe

AST

Automatic synchronization technology



LOW-NOISE BLOWER



APRT


Automatic Pressure Regulation Technology



Comfortable

COMF

Pressure relief at the end of exhalation




ONE BUTTON CONTROL



PREHEAT

5 minutes from 23°C to 45 °C



User-friendly

Alarm Function

10 abnormal alarms including low tide volume



SEPARABLE DESIGN



USB STORAGE



SHIMA RAMPERSAD
Managing Director



084 601 6699
072 028 0307
sales@lotusmedicaire.co.za
www.lotusmedicaire